

# SILVER PRODUCTS GUIDELINE

## Hydrogel

**Minimal to No drainage - Stage III or IV - Change up to once per day.**

**HYDROGEL Ag** by **GENTELL®**

- Aloe Vera based hydrating wound gel, that protects the wound bed and enhances the moist environment.

**SilvaSorb®** by **Medline**

- Combines controlled-release, antimicrobial silver with fluid management.
- Reduces bioburden.
- Constant antimicrobial protection - Controlled-release.

## Collagen

**Change daily as needed.**

**PRISMA**

by **Systagenix**  
WOUND MANAGEMENT

- 55% Collagen - 44% ORC - 1% Silver - ORC

**PURACOL® PLUS AG+** by **MEDLINE**

- 100% pure native collagen

## Alginate

**Moderate to Heavy drainage - Stage III or IV - Change up to once per day.**

**ALGICELL™ Ag:**

by **DERMA SCIENCES**

- Provides superior bactericidal action that goes to work almost immediately, with kill rates that approximate 99.99%
- ALGICELL™ Ag, comprised of a proprietary blend of mannuronic (high gelling) and guluronic (high wet-strength) alginic acids, is designed to gel well while keeping its physical integrity for intact removal.

**Maxorb® Extra Ag** by **MEDLINE**

- Highly absorbent, non woven pad composed of a high G calcium alginate, carboxymethylcellulose and ionic silver.
- Maxorb® Extra Ag is biocompatible, non-irritating, non-sensitizing and will not harm new granulation tissue.

**RESTORE**

by **hollisterwoundcare**

- Odor reduction from antibacterial effect.
- Wound fluid absorbs the alginate and forms a gel assisting in maintaining a moist environment and allowing intact removal.
- Effective antimicrobial agent against a broad spectrum of micro-organisms.

**SeaSorb®-Ag** by **Coloplast**

- Absorbency of SeaSorb-Ag is 33% higher than Aquacel AG.
- Provides an effective, sustained release of silver for up to four days.

**SILVERCEL**

by **Systagenix**  
WOUND MANAGEMENT

- Sustained release of silver ions combats bacterial and fungal growth and help reduce infection in the wound.
- Effectively inhibited the growth in vitro of 150 clinically isolated microorganisms, including antibiotic resistant strains.

## Foam

**Moderate to Heavy drainage - Stage III or IV - Change up to 3 times per week.**

**OptiFoam® AG** by **MEDLINE**

- Controlled-release silver technology kills a broad spectrum of microbes without harming healthy tissue for up to 7 days.
- Highly absorbent foam pad has a smooth, micropore surface to help prevent adherence to the wound bed and minimize trauma during dressing changes.

**PolyMem® Silver™** by **FERRIS**

- Continuous silver reservoir allows effective antimicrobial release. Only silver dressing that contains a cleanser, moisturizer, a super absorbent starch copolymer, and nanocrystalline silver particles.

1-800-604-6068

